

COMPRESSOR NEBULIZER

WOOFY THE DOG

LVA3085DOG

vive
vivehealth.com



OWNER'S MANUAL

TABLE OF CONTENTS

What's Included	3
Intended Use	3
Environment	4
User Profile	4
Contraindications	4
Side Effects	4
About The Compressor Nebulizer	5
What Is Compressor Nebulizer?	5
How Does The Nebulizer Kit Work?	5
Important Use & Safety Information	5
Operating Instructions	6
Cleaning & Maintenance	8
Daily Cleaning	8
Replacement Schedule	9
Storage Tips	9
Important Safety Notes	9
Disposal	9
Troubleshooting	10
Technical Specifications	11
Electromagnetic Compatibility	12
Important Safety Warnings	12
Medical Safety	12
Device Safety	12
Usage Requirements	13
Maintenance & Care	13
Emergency Procedures	13
Signs & Symbols	14
Warranty Information	15
Manufacturer Information	15
Important Safety Notice	16

COMPRESSOR NEBULIZER

WHAT'S INCLUDED



- Compressor Nebulizer Unit (Dog)
- Nebulizer Cup
- Air Tube
- Mouthpiece
- Child Mask
- Replacement Air Filters (x5)
- User Manual

INTENDED USE

The compressor nebulizer is designed to convert liquid medication into a fine aerosol mist for effective delivery to the respiratory system. This device is specifically intended for home respiratory treatment under adult supervision, as prescribed by a healthcare provider. Ideal for delivering medication to treat asthma, allergies, bronchitis, and other respiratory conditions in patients.

ENVIRONMENT

The device is suitable for inhalation at home. Medicines should only be inhaled following the doctor's advice. Besides, inhalation should be performed in a calm and relaxed environment, which should be slow and deep to enable the medicine to reach the small bronchi deep in the lungs. Exhalation should be normal.

USER PROFILE

Caregiver or patient under the guidance of qualified medical experts for home treatment.

The user should understand the general operation of the nebulizer and the content of the instruction manual.

CONTRAINDICATIONS

This device is only for use by awake patients in a home environment. It should not be used during sleep or on unconscious patients. The device is not intended for life support or emergency situations. Consult with a healthcare provider before use if the child has any severe respiratory conditions, recent respiratory surgery, or specific medication requirements.

SIDE EFFECTS

- Some sensitive patients may feel chest tightness and suffer anoxia during aerosol treatment.
- Respiratory or digestive tract infections.

- Hypoxia and discomfort
- Cough
- Circumoral paleness
- Tachycardia

ABOUT THE COMPRESSOR NEBULIZER

WHAT IS A COMPRESSOR NEBULIZER?

A compressor nebulizer is a medical device that converts liquid medication into a fine mist (aerosol) for inhalation. Using the Venturi principle, it combines compressed air with liquid medication to create droplets suitable for effective delivery to the respiratory system.

HOW DOES THE NEBULIZER KIT WORK?

It generates a stream of air that passes through the nebulizer cup. This air mixes with the liquid medication, creating a fine mist that the patient inhales through the mouthpiece or mask. The process ensures optimal medication delivery to the lungs and is especially helpful for kids.

IMPORTANT USE & SAFETY INFORMATION

- Clean and disinfect the nebulizer kit components completely before first use!
- Before each treatment check that:
 - All connections are secure.
 - Tube and medicine cup are tightly connected.
 - Device works by testing without medicine.
- Replace the air filter monthly.

- Replace the nebulizer kit mask if damaged or after 6 months.
- Use a new nebulizer kit when switching medicines.
- Only use prescribed nebulizer medications.
- Never use essential oils, cough medicines, gargles, skin drops, or bath solutions, these can damage your device and cause breathing problems.
- Always check with your doctor about which medications are safe to use.

OPERATING INSTRUCTIONS

1. Open the nebulizer cup by twisting counter-clockwise.

REUSABLE



DISPOSABLE



2. Add the prescribed medication (max. capacity: 8mL).

REUSABLE



DISPOSABLE



NOTE: If the medication volume is less than 2mL, dilute it with saline solution to a minimum of 4mL or as directed by a doctor.

Use medication on the advice of a physician only, and check the appropriate inhalation duration and quantity for your own needs.

3. Close the lid by twisting clockwise.

REUSABLE



DISPOSABLE



4. Choose your inhalation device.



a. **Mouthpiece:** Best option if your child can use it.



b. **Mask:** Good option for younger children and is hand's free.

5. Keep the nebulizer cup upright at all times and have your child breathe in slowly and deeply.
6. Treatment will end once the medicine has evaporated completely, disconnect the air tubes and clean all components (refer to “Cleaning & Maintenance” Section)

CLEANING & MAINTENANCE

DAILY CLEANING

After each use:

Reusable nebulizer kit can be cleaned daily in dishwasher on the top shelf.

- Disassemble mouthpiece or mask from cup, open nebulizer by turning cap counterclockwise and remove the top.

Weekly Disinfection (Choose one):

- Soak in 1:3 vinegar-water solution for 30 minutes, rinse, dry.
- Boil heat-resistant parts for 5 minutes, cool, dry.

Wipe main unit with slightly damp cloth only, never submerge

IMPORTANT: Ensure all parts are completely dry before storage.

REPLACEMENT SCHEDULE

Part	When to Replace
Air Filters	Every month <ul style="list-style-type: none"> • Locate air filter cover and remove it • Pinch the air filter to grab and remove it:
Nebulizer Kit*	<ul style="list-style-type: none"> • Every 3 months with heavy use • Every 6 months with light use • Immediately if damaged

*Nebulizer kit includes mask, air tubes, and medicine cup.

STORAGE TIPS

- Keep in dry place away from sunlight.
- Store in original packaging when possible.
- Avoid bathroom or damp areas.
- Don't store with wet items.

IMPORTANT SAFETY NOTES

- Never use dishwasher.
- Keep unit unplugged during cleaning.
- Never spray liquid into vents.
- Check with doctor if you have special medical conditions (like cystic fibrosis).

DISPOSAL

- Don't dispose with household waste.
- Follow local electronics recycling guidelines.
- Contact local authorities for proper disposal methods.

TROUBLESHOOTING

Problem/Question	Possible Cause/Remedy
The nebulizer produces little to no aerosol	<ul style="list-style-type: none"> • Too much or too little medicine in the nebulizer cup (Minimum: 0.5mL; Maximum: 8mL) • Check and ensure that the nozzles inside the nebulizer cup are not blocked. If necessary, clean the nozzle (e.g. by rinsing). Use the nebulizer cup again. <p>CAUTION: Only poke through the fine holes from the underside of the nozzle.</p> <ul style="list-style-type: none"> • The nebulizer cup is not held in an upright position. • The medicine solution is unsuitable for nebulizing (e.g. too thick). The medicine solution should be specified by the doctor.
Output is too low	Tube kinked, filter blocked, too much inhalation solution.
Which medicines are suitable for inhalation?	Consult your doctor. Essentially all medicines that are suitable and approved for device-based inhalation can be inhaled.
Some inhalation solution remains in the nebulizer	This is normal, and occurs for technical reasons. Stop inhalation if you hear a notable difference in the noise of the nebulizer.
What special steps should be taken for babies and children?	<p>For babies, use the child mask. The mask should cover the mouth and nose to guarantee effective inhalation. For babies, use the child mask. The mask should also cover the mouth and nose.</p> <p>Nebulizing is not appropriate for sleeping people because there is not enough medicine reaching the lungs.</p> <p>Note: Children and babies should only use the device under the help and supervision of an adult. Never leave a child alone with the nebulizer.</p>
Inhalation with the mask takes a longer time	This is for technical reasons. Less medicine is inhaled per breath through the mask holes than using the mouthpiece. The aerosol is mixed with ambient air through the holes in the mask.
Does each person need their own nebulizer?	This is absolutely necessary for hygiene reasons.

TECHNICAL SPECIFICATIONS

Name	Compressor Nebulizer
Model No.	Dog: MCN-S600MB
Power	AC 230V, 50HZ
Power Consumption	below 65W
Nebulizer Flow	≥5L/min
Compressor Flow	≥8L/min
Medicine Container Capacity	8ml
Particle Size	≥60% within 0.5 ~5um range
MMAD	2.80±25%µm
Average Nebulization Rate	0.25 ml/min
Residual Volume	below 1 ml
Noise Level	below 58dB(A)
Max Pressure	≥ 29 Psi / ≥ 200 Kpa
Operation Pressure Range	≥ 12 Psi / ≥ 83 Kpa
Operating Temperature Range	10°C to 40°C (50°F to 104°F)
Operation Humidity	10% to 95%RH
Storage Temperature Range	-20°C to 70°C (-4°F to 158°F)
Storage Humidity Range	10% to 95%RH
Product Weight	2.87lbs
Dimension (LxWxH)	4.72 x 7.87 * 8.47in
Standard accessories	Filters (5), Air Tube, Medicine Container (8 mL), Mouthpiece, Child Mask

ELECTROMAGNETIC COMPATIBILITY

- The device complies with current specifications about to electromagnetic compatibility, and is suitable for use in all premises, including those designated for private residential purposes.
- The radio frequency emissions of the device are extremely low and, in all probability, do not cause any interference with other devices in the proximity.
- It is recommended that you do not place the device on top of or close to other electrical devices, move the unit or connect it to a different socket.
- Radio equipment may affect the operation of this device.

IMPORTANT SAFETY WARNINGS

MEDICAL SAFETY

- Only use medicines prescribed by your doctor.
- Follow your doctor's instructions for dosage and frequency.
- Stop use immediately if you experience discomfort or breathing problems.
- Not for life support or emergency medical treatment.
- Keep medication and mist away from eyes.

DEVICE SAFETY

- Inspect device before each use for damage.
- Only use with approved nebulizer accessories.
- Keep device dry and away from water.
- Unplug immediately if device makes unusual sounds.

- Don't use near flammable gases or high oxygen concentrations.
- Keep away from heat sources.
- Don't block ventilation openings.
- Use correct voltage as marked on device.

USAGE REQUIREMENTS

- Adult supervision required for children and disabled persons.
- Children under 5 must use mask instead of mouthpiece.
- Each person needs their own nebulizer kit for hygiene.
- For household use only, not commercial/clinical use.
- Use only in well-ventilated areas.

MAINTENANCE & CARE

- Disconnect power after each use.
- Never repair unit yourself - contact customer service.
- Keep unit clean and dust-free.
- Store in specified conditions away from moisture.

EMERGENCY PROCEDURES

If device overheats:

1. Turn off immediately.
2. Unplug from outlet.
3. Contact customer service if unusual noise persists.

For technical support or replacement parts, contact Vive Health customer service.

SIGNS AND SYMBOLS

	Warning/Danger: Improper use might cause serious injuries
	Following operating instructions
	Manufacturer information: see the coverage page
	Serial number
IP21	IP code of the device: the grade of this device against ingress of solid foreign objects - 12.5mm diameter (and against access to hazardous parts with finger); the grade of waterproof is protected against vertically falling water drops.
 1639	Complies with the European Medical Device Directive (93/42/EEC and amended Directive 2007/47/EC. Notified Body is SGS Belgium NV.
	Disposal in accordance with Directive 2002/96/EC (WEEE)
	The notices should always be followed
	Applied part of type BF
	Type of protection against electric shock: Class II equipment
	Batch code

WARRANTY INFORMATION

Vive Health warrants this product to be free from manufacturing and material defects for a period of 5 years from the date of purchase under normal use. This warranty extends only to the original purchaser and does not cover:

- Damage from improper use, abuse, or modification.
- Normal wear and tear.
- Battery or replaceable components.
- Products purchased from unauthorized sellers.
- Commercial, institutional, or rental use.

To obtain warranty service, contact Vive Health customer service with proof of purchase and a description of the defect. Vive Health will repair or replace defective products at our discretion.

MANUFACTURER INFORMATION



Vive Health

8955 Fontana Del Sol Way
Naples FL 34109

service@vivehealth.com
1-800-487-3808
vivehealth.com

IMPORTANT SAFETY NOTICE!!

If you experience any serious incident related to the use of this device, immediately discontinue use and report it to:

1. The manufacturer: Vive Health (8955 Fontana Del Sol Way, Naples FL 34109)
2. Contact your national health authority
3. Your healthcare provider.

A serious incident is any incident that directly or indirectly led, might have led, or might lead to:

- Death of a patient, user, or other person
- Temporary or permanent serious deterioration of health
- A serious public health threat

For non-serious concerns, contact Vive Health customer service:

service@vivehealth.com

1-800-487-3808

vivehealth.com



HAVE MORE QUESTIONS?

Feel free to connect with our customer service team by phone, e-mail, or chat on our website

service@vivehealth.com
1-800-487-3808
vivehealth.com

Distributed by



8955 Fontana Del Sol Way
Naples, FL 34109
Made in China