

OWNER'S MANUAL
360 WALKER BOOT TALL
SKU: SUP3091BLKT



SUP3091BLKT

360 WALKER BOOT TALL

PUTTING THE WALKER ON



For a video demonstration,
check out coretechortho.link/1fc

1. Unfasten all straps and remove the soft padded liner.



2. Remove the two hook and loop strap protectors from the inside, pull up and discard. This will allow the liner to secure into place.



3. Place the padded liner back into the boot and insert your foot. Wrap the inner ankle pad over the top of your ankle.



4. Secure the bottom, then the top strap to a comfortable fit.



5. Secure the outer straps from bottom to top and tighten to a comfortable fit. Do NOT overtighten.



SUP3091BLKT

6. If the boot is uncomfortable, loosen the straps and use the extra pads included to pad the inner liner or straps.



7. Press the valve to inflate to your desired compression.



8. Press the metal button above the pump to deflate the bladder.



 **coretech**®
by **vive health**

Distributed by

**vive
health**

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com