

WALKER FOREARM PLATFORM


QUICK START GUIDE



To see all FAQs in one place, visit vhealth.link/qje

WHAT'S INCLUDED

- (A) 1x Upright Forearm Base
- (B) 2x Attachment Brackets
- (C) Hand Grip Extension



For a video demonstration, visit vhealth.link/vg8

1. Slide both of the Attachment Brackets (B) onto the post of the Forearm Base (A)



2. Attach the brackets onto the upper and lower side posts of your walker. You may have to loosen the knobs to get them to fit onto the posts. Make sure when doing this that the Forearm Bases armrest is hanging over the walker and not facing outward. This is very important as if you incorrectly face the armrest outward it will be dangerous and much easier to tip.

(**WARNING:** If your walker does not have an upper or lower post on the side of your walker to safely attach both brackets then this post was not made if your walker and you do NOT own the STANDARD type walker this item was made for.)



vive®

Distributed by

vive
health

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com
Made in China

3. After fitting the brackets onto both the upper and lower posts of your walker, find a height most suitable for the user and then tighten the brackets down using all four knobs on the two brackets.



4. Insert Hand Grip Extension (C) into the front of the Forearm Base and tighten down the two knobs underneath after you find a length and position most suitable for the user.



5. When Finished the unit should look like this:



WARNINGS ⚠️

- Max user weight should not exceed 300lbs.
- Do not use if this item doesn't properly fit your style of walker. This item will fit most traditional style walkers.
- Carefully follow all instructions and make sure your Walker Forearm Platform looks exactly the way it is set up in step 5 with the forearm section draped inward to the inside of the walker.
- Make sure your arm is strapped securely into the forearm rest before use