
ROHO® DRY FLOATATION® Cushions Operation Manual

Single Compartment®, Sensor Ready®, ENHANCER®, Quadtro SELECT®, Contour SELECT®

permobil

s h a p e f i t t i n g t e c h n o l o g y®

Contact Information

ROHO, Inc.
1501 South 74th St.
Belleville, IL 62223 U.S.A.

U.S.A.:
1-800-736-0925
Fax 1-888-551-3449

Customer Support:
orders.roho@permobil.com

Outside the U.S.A.:
+1 618-277-9150
Fax +1 618-277-6518
intlorders.roho@permobil.com

permobil.com

Warning

Indicates that not following the specified procedure could lead to potentially hazardous conditions resulting in serious injury.

Caution

Indicates that not following the specified procedure could lead to potentially hazardous conditions resulting in minor to moderate injury or damage to the equipment or other property.

Incident Reporting

In the event of a serious incident or injury, contact Customer Support and the competent authority of the State in which the user and/or patient is established.

Contents

Intended Use			2-3
Important Safety Information			
Parts Detail			4
Product Specifications			
Cushion Setup			4-6
Troubleshooting			7
Cleaning and Disinfecting			7-8
Storage, Transport, Disposal and Recycling			8
Warranty			
Español	9	Eesti	109
Français	17	Ελληνικά	117
Deutsch	29	magyar	125
Italiano	37	Íslenska	133
Nederlands	45	Lietuviškai	141
Dansk	53	Polski	149
Svenska	61	Română	157
Norsk	69	Slovenčina	165
Suomi	77	Slovenščina	173
Português	85	Türkçe	181
hrvatski	93	Српски	189
Čeština	101		

Symbols Glossary

	Manufacturer		Fits		Size (width x depth x height)		Front of product
	Manufacturer's catalog number (in UDI - unique device identification)		Consult instructions for use		No weight limit, if sized to user		Back of product
	Model number for customer reference only		Serial number		Medical Device		Bottom of product
	Europe Authorized Representative		Quantity		Date of manufacture		Batch code
	Swiss Authorized Representative		Unique Device Identifier		DO NOT use the ISO FLO Memory Control as a handle; the product may tear.		

Please Note

Supplier: This manual must be given to the user of this product.
Operator (Individual or Caregiver): Before using this product, read the instructions and save for future reference.

Products may be covered by one or more U.S. and foreign patents and trademarks, including CONTOUR SELECT®, DRY FLOATATION®, ENHANCER®, HIGH PROFILE®, ISO FLO®, LOW PROFILE®, MID PROFILE®, QUADTRO SELECT®, ROHO®, Sensor Ready®, shape fitting technology®, Smart Check® and Yellow Rope®. Zytel® is a registered trademark of DuPont.

Product information is changed as needed; current product information is available at permobil.com. Contact Customer Support for previous versions of product information.

ROHO, Inc. is part of Permobil.

© 2007, 2026 Permobil

Rev. 2026-01-09

The following products are featured in this operation manual:

ROHO LOW PROFILE Single Compartment Cushion

ROHO HIGH PROFILE Single Compartment Cushion

ROHO HIGH PROFILE Single Compartment Cushion with Sensor Ready Technology

ROHO ENHANCER Cushion

ROHO CONTOUR SELECT Cushion

ROHO QUADTRO SELECT LOW PROFILE Cushion

ROHO QUADTRO SELECT MID PROFILE Cushion

ROHO QUADTRO SELECT HIGH PROFILE Cushion

Intended Use

Each ROHO DRY FLOATATION product in this manual is an adjustable, air-filled, cellular-design wheelchair support surface. There is no weight limit, yet the cushion must be properly sized to the individual.

Each of the following products utilizes DRY FLOATATION Technology and is intended to conform to an individual's seated shape to provide skin/soft tissue protection, positioning, and an environment to facilitate wound healing.

SELECT and ENHANCER Cushions:

ROHO QUADTRO SELECT LOW PROFILE Cushion (QUADTRO SELECT LOW PROFILE)

ROHO QUADTRO SELECT MID PROFILE Cushion (QUADTRO SELECT MID PROFILE)

ROHO QUADTRO SELECT HIGH PROFILE Cushion (QUADTRO SELECT HIGH PROFILE)

ROHO CONTOUR SELECT Cushion (CONTOUR SELECT)

ROHO ENHANCER Cushion (ENHANCER)

Each of the following products utilizes DRY FLOATATION Technology and is intended to conform to an individual's seated shape to provide skin/soft tissue protection and an environment to facilitate wound healing.

Single Compartment Cushions:

ROHO LOW PROFILE Single Compartment Cushion (LOW PROFILE Single Compartment Cushion)

ROHO HIGH PROFILE Single Compartment Cushion (HIGH PROFILE Single Compartment Cushion)

ROHO HIGH PROFILE Single Compartment Cushion with Sensor Ready Technology (HIGH PROFILE Sensor Ready Cushion)

The Sensor Ready Cushion is compatible with the Smart Check Cushion Monitoring Device.

Contraindication: Based on clinical, scientific or engineering evidence, these products may not be suitable for individuals who require accommodation of pelvic asymmetry greater than is specified below, per product:

1 in. (2.5cm) - QUADTRO SELECT LOW PROFILE, ENHANCER, CONTOUR SELECT, LOW PROFILE Single Compartment Cushion

2 in. (5cm) - QUADTRO SELECT MID PROFILE

3 in. (7.5cm) - QUADTRO SELECT HIGH PROFILE, HIGH PROFILE Sensor Ready Cushion, HIGH PROFILE Single Compartment Cushion

A clinician who is experienced in seating and positioning should determine whether the product is appropriate for the individual's particular seating needs. The product's intended use statement alone should not be relied upon to make that determination.

The medical products featured in this document are intended to be only one part of an overall care regimen that includes all seating and mobility equipment and therapeutic interventions. A clinician who is experienced in seating and positioning should identify this care regimen after evaluating the individual's physical needs and overall medical condition.

A clinician should also evaluate for visual, reading, and cognitive impairment to determine the need for caregiver assistance or the need for other assistive technology, such as large-print instructions, to ensure proper product use.

Important Safety Information

Warnings:

- Skin/soft tissue breakdown can occur due to a number of factors, which vary by individual. Check skin frequently, at least once a day. Redness, bruising, or darker areas (when compared to normal skin) may indicate superficial or deep tissue injury and should be addressed. If there is any discoloration to skin/soft tissue, STOP USE immediately. If the discoloration does not disappear within 30 minutes after disuse, immediately consult a healthcare professional.
- DO NOT place obstructions between the individual and the cushion. The cushion and the cover MUST be compatible sizes and MUST be used as directed in this manual. Except for the compatible covers listed in "Product Specifications" in this manual, placement of any items between the individual and the cushion: 1) may reduce or eliminate the benefits of the cushion, increasing the risk to skin or to other soft tissue, and 2) may cause the individual to become unstable and vulnerable to falling.
- DO NOT use the cushion on top of, or in conjunction with, any other products or materials, except as indicated in this manual. Using non-compatible products or materials, whether inside or outside the cover: 1) may reduce or eliminate the benefits of the cushion, increasing the risk to skin or to other soft tissue, and 2) may cause the individual to become unstable and vulnerable to falling.

Warnings, continued:

- If unable to perform any tasks described in this manual, seek assistance by contacting the clinician, equipment provider, distributor, or Customer Support.
- DO NOT expose the cushion to high heat, open flames or hot ashes. Testing or certification claims, including for flammability, may no longer apply to this device when it is combined with other products or materials. Check testing and certification claims for all products used in combination with this device.
- The cushion surface will acclimate to the surrounding temperature. Take appropriate precautions, especially when the cushion will contact unprotected skin.
- Use caution when handling the cushion: Metal grommets on the corners may be a finger trap hazard.
- DO NOT use the cushion as a water flotation device (e.g., a Life Preserver). It will NOT support the body in water.

Important Safety Information, continued**⚠ Cautions:**

- Keep the cushion away from sharp objects.
- Changes in altitude may require adjustment to the cushion.
- If the cushion has been in temperatures lower than 32°F (0°C) and exhibits unusual stiffness, allow the cushion to warm to room temperature.
- DO NOT use a pump, cover, or repair kit other than those compatible with the product.
- DO NOT modify the product or use non-compatible accessories or components, which may result in product damage or failure and will void the warranty.
- DO NOT allow the cushion to come into contact with oil-based lotions or lanolin. They may degrade the material.
- Prolonged exposure to ozone may degrade materials used in the cushion, affect the performance of the cushion, and void the product warranty.
- Periodically, check the cover and components for damage and replace if necessary. See "Troubleshooting".
- DO NOT use any of the following as a handle for carrying or pulling the cushion: inflation valve(s), hose(s), quick disconnect (Sensor Ready Cushion), or the ISOFLO Memory Control (SELECT Cushion). Carry the cushion by the cushion base or by the carrying handle on the Standard/SELECT Cover. Using any other components to carry or pull the cushion may damage the product.

⚠ Warnings - Motor Vehicle Transportation:

- Failure to pay attention to these warnings could result in severe injury to the individual in the wheelchair or to others.
- Whenever possible, transfer out of the wheelchair and into a manufacturer-installed vehicle seat without the cushion, and use the vehicle's crash-tested occupant restraint system.
- The cushion must be correctly and securely installed, used, and maintained according to the instructions and warnings provided by the cushion manufacturer.
- When the wheelchair is used as the vehicle seat, the cushion should only be used with wheelchair frames that are equipped with securement points specified in ANSI/RESNA WC-4:2017, section 19 (WC19) or ISO 7176-19, and that have been successfully crash tested on a wheelchair in the forward-facing position.
- The wheelchair with the cushion must be used with an effective wheelchair securement system and a properly positioned, crash-tested pelvic and shoulder-belt restraint, or Wheelchair Tiedown and Occupant Restraint System (WTORS), following the manufacturer's instructions.
- If the cushion has been involved in an accident during transport, check the cushion for damage and replace if necessary. See the "Troubleshooting" section in this manual.
- For more information about safely transporting wheelchair users, refer to: <https://wc-transportation-safety.umtri.umich.edu/ride-safebrochure>.
- A ROHO QUADTRO SELECT Cushion was included in a wheelchair seating system that was dynamically tested by a third-party testing facility for use in a motor vehicle. The system tested meets all applicable criteria for wheelchair seating systems proposed in Section 5.1 of ANSI/RESNA WC-4:2017 Section 20 "Wheelchair Seating Systems for Use in Motor Vehicles."
- Test results apply to seating systems intended for use by adults and children with a total body mass greater than 23 kg (51 lb, or approximately 6 years of age) who are transported facing forward in their wheelchairs in all types of motor vehicles.
- To determine whether motor vehicle third-party testing has been performed with a specific seating system combination, refer to all provided equipment operation manuals or contact the equipment manufacturer.


Parts Detail

Package contents: cushion, cover*; repair kit, hand pump, operation manual.

* Instructions are provided with the cover.

To order replacement parts, contact the product provider or Customer Support.

The serial number label is located on the bottom of the cushion.

Product Specifications

Cushion Dimensions and Weights

	SINGLE COMPARTMENT Cushions		QUADTRO SELECT Cushions			Contoured Cushions	
	LOW PROFILE	HIGH PROFILE**	LOW PROFILE	MID PROFILE	HIGH PROFILE	CONTOUR SELECT	ENHANCER
Height	2¼ in. (5.5cm)	4¼ in. (10.5cm)	2¼ in. (5.5cm)	3¼ in. (8.5cm)	4¼ in. (10.5cm)	2¼ to 4¼ in. (5.5 to 10.5cm)	2¼ to 4¼ in. (5.5 to 10.5cm)
Average Weight*	2.5 lbs. (1.1 kg)	3.25 lb. (1.5 kg)	2.5 lbs. (1.1 kg)	3.0 lbs. (1.4 kg)	3.25 lb. (1.5 kg)	3.5 lbs. (1.6 kg)	3.5 lbs. (1.6 kg)
Min. Width	6 cells; 11½ in. (29cm)	6 cells; 11½ in. (29cm)	5 cells; 11 in. (28cm)	7 cells; 14.5 in. (36.5cm)	6 cells; 12¾ in. (32.5cm)	7 cells; 14½ in. (36.5cm)	8 cells; 14¾ in. (37.5cm)
Min. Depth	6 cells; 11½ in. (29cm)	6 cells; 11½ in. (29cm)	6 cells; 11¼ in. (30cm)	8 cells; 15¼ in. (38.5cm)	6 cells; 11½ in. (30cm)	8 cells; 15¼ in. (38.5cm)	8 cells; 14¾ in. (37.5cm)
Max. Width	13 cells; 23½ in. (59.5cm)	15 cells; 26¼ in. (68cm)	13 cells; 24¼ in. (62.5cm)	11 cells; 21¼ in. (54cm)	14 cells; 26¼ in. (67cm)	13 cells; 24½ in. (62.5cm)	11 cells; 20 in. (50.5cm)
Max. Depth	12 cells; 21¾ in. (55cm)	14 cells; 23½ in. (59.5cm)	12 cells; 22 in. (56cm)	12 cells; 22 in. (56cm)	14 cells; 23¾ in. (60cm)	12 cells; 22 in. (56cm)	12 cells; 21¼ in. (55cm)

* based on 9 x 9 cell cushions: 16½ in. x 16½ in. (42cm x 42cm)

1 in. = 2.54 cm = 25.4 mm

** applies to the Sensor Ready HIGH PROFILE Cushion, except Max. Depth is 12 cells; 21¼ in (55cm).

Note: Contact Customer Support for specific cushion sizes.

Covers and Accessories Compatible with ROHO DRY FLOATATION Cushions:

Not available in all markets. Contact Customer Support.

- ROHO Heavy Duty Cushion Cover
- ROHO Standard Cushion Cover
- ROHO SELECT Cushion Cover
- ROHO Planar Solid Seat Insert
- ROHO Contour Base
- Smart Check Cushion Monitoring Device

Cushion materials: cushion: neoprene rubber, made without latex; inflation valve: Zytel; ISO FLO Memory Control: polypropylene, polycarbonate, aluminum; Quick Disconnect: acetal

Cover materials - Refer to the care label on the cover.

Cushion expected lifetime: 5 years

SMART CHECK OWNER: The Sensor Ready Cushion is equipped with Sensor Ready Technology and can be used with Smart Check, which is a cushion feedback system. Smart Check can be used to find the individual's inflation level and help check their cushion. For more information, contact Customer Support.

Note: STOP HERE and refer to the cushion setup instructions provided with Smart Check.

Cushion Setup

Warnings:

- DO NOT use a product that is underinflated or overinflated, because 1) the product benefits will be reduced or eliminated, resulting in an increased risk to skin and other soft tissue, and 2) the individual may become unstable and vulnerable to falling. Carefully follow the instructions for inflation, placement and hand check. If the product does not appear to be holding air, or if unable to inflate or deflate the product, see "Troubleshooting". Immediately contact the equipment provider, distributor, or Customer Support if the problem persists.
- When possible, use hook and loop fasteners to further secure the cushion to the seating surface. Hook fasteners are available upon request.
- After setting up the cushion the first time, perform a hand check frequently, at least once a day. When the cushion is used by a different individual or in a different wheelchair, repeat cushion setup. Follow the hand check instructions in this manual.

Cushion Setup is continued on the next page.


Cushion Setup, continued

Inflate, cover, and transfer.

IMPORTANT - For SELECT Cushions:

Before inflating, unlock the ISO FLO Memory Control so that all of the air cells will inflate. Push the green knob to unlock the ISO FLO Memory Control.

Note: The ISO FLO can be accessed through an opening in the front of the ROHO SELECT Cushion Cover (if used).

- 1 Inflate:** Slide the hand pump nozzle over the inflation valve. Open the inflation valve. Inflate the cushion until all of the air cells feel firm.

Close the inflation valve until it comes to a stop. Remove the hand pump. Make sure the inflation valve is closed. Repeat for all inflation valves.

Note: To open, turn the inflation valve at least one full rotation.


- 2 Cover:** Place the cushion inside the cover. Refer to the instructions provided with the cover.
Note: Keeping a cover on the cushion will help to protect the cushion components.

3 Transfer:

- a. Place the cushion on the wheelchair with the air cells facing up and the non-skid surface of the cover facing down. The symbol showing the seated figure on the cover may be used to help orient cushion. (Fig.1)

Note: Consult the prescriber about alternate positions of inflation valve(s) and ISO FLO Memory Control (for SELECT Cushions).

- b. The individual should transfer to the wheelchair and sit in their usual sitting position. Make sure that the cushion is centered underneath the individual. The inflation valve should be in the front, on the left side of the seated individual. (Fig. 2)

Perform a hand check.

- 4** To perform a hand check, use the width of fingers to check the amount of air in the cushion. When the cushion has the correct amount of air*, the individual is both immersed in the cushion and supported by air:

- a. Slide a hand between the cushion and the individual, and feel for the lowest bony prominences. *Suggestion: To help locate the bony prominences while performing a hand check, lift and lower the leg.*
b. Refer to the Hand Check Guide below to adjust the cushion.

* The recommended distance between the lowest bony prominence and the seat surface is $\frac{1}{2}$ - 1 in. (1.5 - 2.5 cm).

IMPORTANT! Perform cushion setup with all valves. SELECT and ENHANCER Cushions may require additional positioning and adjustment. See additional instructions on the next page.


HAND CHECK GUIDE

UNDERINFLATED!

Not able to move fingers.

Not enough air in the cushion.

Risk: Not supported by air.


OVERINFLATED!

Can move fingers too freely.

Too much air in the cushion.

Risk: Not immersed in the cushion.

CORRECT AMOUNT OF AIR!

Able to slightly move fingers.

Immersed in the cushion and supported by air.

If the cushion is UNDERINFLATED:

Add some air. Close the inflation valve. Perform a hand check again. Repeat until the cushion has the correct amount of air.

If the cushion is OVERINFLATED:

Remove some air. Close the inflation valve. Perform a hand check again. Repeat until the cushion has the correct amount of air.

If the cushion has the CORRECT AMOUNT OF AIR, it is ready to use.

Cushion Setup, continued


ENHANCER

- Follow the instructions on the previous page for Cushion Setup. Explanation of inflation valves:

"A" - inflation valve for outer cells/abductor (air cells shown in white)

"B" - inflation valve for posterior contour section (air cells shown in gray)

- Make sure that both inflation valves are closed when setup is complete.

SELECT Cushion Adjustment for a Fixed or Flexible Pelvic Obliquity or Tilt


⚠ Warnings - SELECT Cushions:

- Consult with the prescriber or a clinician about adjustment of inflation levels in SELECT Cushion compartments.
- After adjusting the SELECT Cushion, lock the ISOFLO as instructed. If the ISOFLO is not locked, the inflation level in each compartment will change.

Fixed Pelvic Obliquity or Tilt:

To adjust the QUADTRO SELECT Cushion or the CONTOUR SELECT Cushion to accommodate a fixed pelvic obliquity or a fixed posterior pelvic tilt:

- Follow the Cushion Setup instructions on the previous page to make sure the cushion has the correct amount of air.
- The individual should sit in their normal position while the air transfers through the ISOFLO Memory Control. When the air transfer is complete, push the red knob toward the locked position.
- Perform a hand check for each compartment of the cushion again, as instructed on the previous page.

GREEN KNOB: Unlock

Flexible Pelvic Obliquity or Tilt:

To adjust the QUADTRO SELECT Cushion to correct a flexible pelvic obliquity or a flexible posterior pelvic tilt:

- Follow the Cushion Setup instructions on the previous page to make sure the cushion has the correct amount of air.
- Make sure the ISOFLO Memory Control is in the unlocked position
- Follow the instructions below.

RED KNOB: Lock

Flexible pelvic obliquity:

- Manually correct to the desired pelvic position. (For example, if the left side of the pelvis is lower than the right side, either have the individual lean toward the right, or place downward pressure on the right hip until the desired position is achieved.)
- Push the red knob toward the locked position.
- Perform a hand check as instructed on the previous page.

Flexible posterior pelvic tilt:

- Check to make sure there is approximately 1" (2.5 cm) of space between the individual's bony prominences and the seating surface. Add more air if needed.
- Lift both of the individual's knees to put more downward pressure on the pelvis and to force air to the front of the cushion. This will create a "well" (a lower, concave area) for the pelvis and will help prevent the pelvis from sliding forward.
- Push the red knob toward the locked position.
- Perform a hand check as instructed on the previous page.

Troubleshooting

For additional assistance, contact the equipment provider, distributor, or Customer Support.

The cushion is not holding air.	Inflate the cushion. Inspect the quick disconnect (Sensor Ready Cushion), the inflation valve(s), and hose(s) for damage. Confirm that the inflation valve(s) are completely closed. Look for holes in the cushion. If very small holes or no holes are visible, follow the instructions in the Repair Kit provided with the product.
Not able to inflate or deflate the cushion.	Confirm that the inflation valve(s) are open. Inspect the inflation valve(s) and hose(s) for damage. When inflating SELECT Cushions: Confirm that the green knob of the ISO FLO Memory Control is pushed toward the center of the ISO FLO in the UNLOCKED position so that all chambers will inflate. For the Sensor Ready Cushion: also check the quick disconnect for damage.
The cushion is uncomfortable or unstable.	- Make sure that the cushion is not underinflated or overinflated (See Cushion Setup). - The cushion cells and the top of the cover should face up. - The cushion should be centered beneath the individual. - The cushion should be properly sized to the individual and the wheelchair. If the cushion still feels uncomfortable or unstable after checking all of the above, immediately consult the prescribing healthcare professional.
The cushion slides on the wheelchair.	Make sure that the non-skid bottom of the cushion cover faces away from the individual. Use directional labels on the cushion and on the cover to orient the cushion. When possible, use hook and loop fasteners to further secure the cushion to the wheelchair. Hook fasteners are available upon request. Refer to the instructions provided with the cover.
A component is damaged or missing.	Contact Customer Support for warranty or replacement.

Cleaning and Disinfecting

Note: For cover care instructions, refer to the manual provided with the cover.

Warnings:

- DO NOT allow water or cleaning solution to enter the cushion. Make sure the inflation valve(s) are closed.
- Sensor Ready Cushion: Make sure the Smart Check is detached from the cushion. Water in a cushion with Sensor Ready (SR) technology can result in damage to the Smart Check and cause a false reading when using Smart Check.
- Cleaning and disinfecting are separate processes. Cleaning must precede disinfection. If the product becomes soiled, or between use by different individuals: clean, disinfect, and check the product for proper functioning.
- Clean the product regularly and when it becomes soiled.
- Follow all manufacturer bleach container and germicidal disinfectant safety guidelines.
- Note about sterilization: High temperatures accelerate aging and will damage the product assembly. The products featured in this manual are not packaged in a sterile condition, nor are they intended or required to be sterilized prior to use. If institution protocol demands sterilization: First, follow the cleaning and disinfecting instructions; then, open the inflation valve(s) and use the lowest possible sterilizing temperature, but not to exceed 70°C (158°F), for the shortest time possible. DO NOT autoclave.

Cautions:

- Use of the following may cause damage to the cushion: abrasives (steel wool, scouring pads); caustic, automatic dishwasher detergents; cleaning products that contain petroleum or organic solvents, including acetone, toluene, Methyl Ethyl Ketone (MEK), naphtha, dry cleaning fluid, adhesive remover; oil-based lotions, lanolin, ozone gas; ultraviolet light cleaning methods.
- DO NOT expose to direct sunlight.

Recommendation: When washing by hand, use room-temperature water.


Note: Deflated air cells may be easier to thoroughly clean. To deflate the air cells: with the inflation valve(s) open, roll up and/or compress the cushion, and then completely close the inflation valve(s).

Hand Wash and Hand Disinfect the Cushion (Preferred Method)

Hand Wash the Cushion: Remove the cover, close the inflation valve(s), and place in a large sink. Hand wash, using mild liquid hand soap, hand dishwashing liquid, laundry detergent, or a multipurpose detergent (follow product label instructions). Gently scrub all cushion surfaces, using a soft plastic bristle brush, a sponge or a cloth. Rinse with clean water. Air dry.

Hand Disinfect the Cushion: Follow the hand wash instructions, using 1 part household liquid bleach per 9 parts water. Keep the product wet with the bleach solution for 10 minutes. Rinse thoroughly with clean water. Air dry.

 **Caution:** Rinse thoroughly. Detergent residue may cause air cells to stick together.

Cleaning and Disinfecting, continued

Machine Wash and Disinfect the Cushion

⚠ IMPORTANT: Machine washing can cause wear and damage, including cuts and tears to the cushion and the detachment of required product information. Washing and disinfecting the cushion by hand and air drying is strongly recommended.

⚠ Warnings:

- DO NOT allow water or cleaning solution to enter the cushion. Make sure the inflation valve(s) are closed.
- Sensor Ready Cushion: Make sure the Smart Check is detached from the cushion. Water in a cushion with Sensor Ready (SR) technology can result in damage to the Smart Check and cause a false reading when using Smart Check.
- High water temperature may cause cushion discoloration and label detachment.
- DO NOT expose cushions to temperatures above 70°C (158°F), which may damage cushion components.
- For machine washing or disinfecting, make sure the washing machine capacity is sufficient. A washing machine that isn't large enough may not remove all soil or contaminants from a cushion.
- Using a damaged cushion will reduce or eliminate the benefits of the cushion and increase risk to skin and other soft tissue.

Prepare the Cushion for Machine Cleaning:

- 1) Remove the cover and close the inflation valve(s). Examine the cushion for damage, and patch any holes (refer to Repair Kit instructions provided with the product).
- 2) Remove stains or odors as needed: Place the cushion in warm water and a disinfectant cleaning agent. Gently scrub, using a soft plastic bristle brush, and rinse with clean water.
- 3) Deflate the cushion: Open the inflation valve and roll up the cushion, starting at the corner opposite the valve. After having removed as much air as possible, close the inflation valve.
- 4) Inspect for any sharp edges in the washing machine that may damage the cushion.

Machine Wash the Cushion: Place the cushion in a laundry bag and then into the washing machine. When folding the cushion, make sure the air cells face outward. Use warm water (85° - 105°F / 30° - 40°C) and laundry detergent or machine- and cushion-safe cleaner/biocide, placed in the detergent dispenser. Gentle cycle: Use a machine setting that does not include a high-speed spin cycle.

Machine Disinfect the Cushion: Follow the machine washing instructions and use hot water 60°C (140°F). Use biocide (refer to the manufacturer instructions). Rinse with clean water.

⚠ Caution: Rinse thoroughly. Detergent residue may cause air cells to stick together.

Drying the Cushion: Remove the cushion from the laundry bag and air dry. For faster drying, use a hand-held hair dryer, hold no closer than 6 in. (15 cm) from the cushion, and use a sweeping motion.

⚠ Warning: Overheating the neoprene rubber can cause permanent product damage.

Care Symbols

Hand wash.	Line dry in the shade.	Bleach (1 part bleach: 9 parts water).	Machine wash, warm, gentle.
Machine wash hot, normal, at the temperature shown.	Use a handheld dryer.		

Storage, Transport, Disposal and Recycling

Storage or Transport: Clean and disinfect the product. Open the inflation valve(s). Deflate the product. Place in a container that protects the product from moisture, contaminants and damage. After storing or transporting, inspect the product to ensure it was not damaged (see "Troubleshooting"), and follow instructions to prepare the product for use.

Disposal: The components of the products in this manual are associated with no known environmental hazards when properly used and when disposed of in accordance with all local/regional regulations. At the end of its useful life, treat as healthcare waste and dispose of according to local healthcare-waste guidelines. Incinerate with other healthcare waste. Incineration must be performed by a qualified, licensed waste management facility.

Recycling: Contact a local recycling agency to determine recycling options for the product.

Warranty

This product is covered by a warranty against manufacturer defects (SELECT Cushions: 36 months; all other DRY FLOATATION products: 24 months) with the exception of product covers that are eligible for warranty up to 6 months.

If a manufacturing defect should occur, discontinue use immediately. Wear from regular use, misuse, damage by an accident or damage from "an act of Nature" (e.g., flood, tornado, earthquake and/or fire) is not considered a manufacturing defect.

For warranty details, please visit permobil.com or contact Customer Support to request a printed copy of the full warranty policy.

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Cuscino ROHO LOW PROFILE Single Compartment

Cuscino ROHO HIGH PROFILE Single Compartment

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Cuscini Single Compartment

Cuscino ROHO LOW PROFILE Single Compartment (cuscino LOW PROFILE Single Compartment)

Cuscino ROHO HIGH PROFILE Single Compartment (cuscino HIGH PROFILE Single Compartment)

Cuscino ROHO HIGH PROFILE Single Compartment con tecnologia Sensor Ready (cuscino HIGH PROFILE Sensor Ready);

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1 pollice (2,5 cm) - QUADTRO SELECT LOW PROFILE, ENHANCER, CONTOUR SELECT, LOW PROFILE Single Compartment

2 pollici (5 cm) - QUADTRO SELECT MID PROFILE

3 pollici (7,5 cm) - QUADTRO SELECT HIGH PROFILE, HIGH PROFILE Sensor Ready, HIGH PROFILE Single Compartment Cushion

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ROHO LOW PROFILE Single Compartment-dynan

ROHO HIGH PROFILE Single Compartment-dynan

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Single Compartment-dynor:

ROHO LOW PROFILE Single Compartment-dynan (LOW PROFILE Single Compartment-dynan)

ROHO HIGH PROFILE Single Compartment-dynan (HIGH PROFILE Single Compartment-dynan)

ROHO HIGH PROFILE Single Compartment-dynan med Sensor Ready-teknologi (HIGH PROFILE Sensor Ready-dynan)

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2,5 cm (1 tum) – QUADTRO SELECT LOW PROFILE, ENHANCER, CONTOUR SELECT, LOW PROFILE Single Compartment-dyna

5 cm (2 tum) – QUADTRO SELECT MID PROFILE

7,5 cm (3 tum) – QUADTRO SELECT HIGH PROFILE, HIGH PROFILE Sensor Ready-dyna, HIGH PROFILE Single Compartment-dyna

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Jastuci Single Compartment:

Jastuk ROHO LOW PROFILE Single Compartment (jastuk LOW PROFILE Single Compartment)

Jastuk ROHO HIGH PROFILE Single Compartment (jastuk HIGH PROFILE Single Compartment).

Jastuk ROHO HIGH PROFILE Single Compartment s tehnologijom Sensor Ready (jastuk HIGH PROFILE Sensor Ready).

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2,5 cm (1 in) – QUADTRO SELECT LOW PROFILE, ENHANCER, CONTOUR SELECT, jastuk LOW PROFILE Single Compartment

5 cm (2 in) – QUADTRO SELECT MID PROFILE

7,5 cm (3 in) – QUADTRO SELECT HIGH PROFILE, jastuk HIGH PROFILE Sensor Ready, jastuk HIGH PROFILE Single Compartment

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Poduszka ROHO LOW PROFILE Single Compartment

Poduszka ROHO HIGH PROFILE Single Compartment

Poduszka ROHO HIGH PROFILE Single Compartment z technologią Sensor Ready

Poduszka ROHO ENHANCER

Poduszka ROHO CONTOUR SELECT

Poduszka ROHO QUADRO SELECT LOW PROFILE

Poduszka ROHO QUADRO SELECT MID PROFILE

Poduszka ROHO QUADRO SELECT HIGH PROFILE

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Poduszki jednosegmentowe (Single Compartment):

Poduszka ROHO LOW PROFILE Single Compartment (LOW PROFILE Single Compartment)

Poduszka ROHO HIGH PROFILE Single Compartment (HIGH PROFILE Single Compartment)

Poduszka ROHO HIGH PROFILE Single Compartment z technologią Sensor Ready (HIGH PROFILE Sensor Ready)

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2,5 cm (1 cal) — poduszki QUADTRO SELECT LOW PROFILE, ENHANCER, CONTOUR SELECT, LOW PROFILE Single Compartment

5 cm (2 cale) — poduszka QUADTRO SELECT MID PROFILE

7,5 cm (3 cale) — poduszki QUADTRO SELECT HIGH PROFILE, HIGH PROFILE Sensor Ready, HIGH PROFILE Single Compartment

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Јастуци Single Compartment:

Јастук ROHO LOW PROFILE Single Compartment (јастук LOW PROFILE Single Compartment)

Јастук ROHO HIGH PROFILE Single Compartment (јастук HIGH PROFILE Single Compartment)

Јастук ROHO HIGH PROFILE Single Compartment са технологијом Sensor Ready (јастук HIGH PROFILE Sensor Ready)

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2,5 cm (1 in) – јастук QUADTRO SELECT LOW PROFILE, ENHANCER, CONTOUR SELECT, LOW PROFILE Single Compartment

5 cm (2 in) – QUADTRO SELECT MID PROFILE

7,5 cm (3 in) – јастук QUADTRO SELECT HIGH PROFILE, HIGH PROFILE Sensor Ready, HIGH PROFILE Single Compartment

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Permobil EMEA AB
Permobil LTD