



# WHEELCHAIR SEATING SIZE & WEIGHT CHART

2520 Pegasus Drive  
Bakersfield, CA 93308

Ph. (661) 393-4877

[www.kernmedicalsupply.com](http://www.kernmedicalsupply.com)

	Up Thru 5'2" Hemi Height				1. Patient weight and height, look up seat width							
	< 5'0"	5'0"	5'2"	5'4"	5'6"	5'8"	5'10"	6'0"	6'2"	6'4"	6'6"	> 6'6"
<100lbs	14"/16"	14"/16"	14"/16"	16"	16"	16"	16"	16"	16"	16"	16"	16"
100lbs	14"/16"	14"/16"	16"	16"	16"/18"	16"	16"	16"	16"	16"	16"	16"
120lbs	18"	18"	18"	18"	16"/18"	16"/18"	16"	16"	16"	16"	16"	16"
140lbs	18"/20"	18"	18"	18"	18"	18"	18"	18"	18"	18"	16"	16"
160lbs	20"	20"	18"/20"	18"/20"	18"	18"	18"	18"	18"	18"	18"	18"
180lbs	20"	20"	20"	20"	20"	18"/20"	18"/20"	18"	18"	18"	18"	18"
200lbs	22"	22"	20"	20"	20"	20"	20"	18"/20"	18"/20"	18"/20"	18"	18"
220lbs	22"	22"	22"	20"	20"	20"	20"	20"	20"	20"	20"	20"
240lbs	24"	24"	22"/24"	22"/24"	20"	20"	22"	22"	22"	22"	22"	22"
260lbs	24"	24"	24"	24"	22"/24"	22"/24"	22"	22"	22"	22"	22"	22"
280lbs	26"	26"	26"	24"	24"	24"	22"/24"	22"	22"	22"	22"	22"
300lbs	26"	26"	26"	26"	26"	24"	24"	24"	24"	22"	22"	22"
320lbs	28"	28"	28"	26"/28"	26"	26"	24"	24"	24"	24"	24"	24"
340lbs	28"	28"/30"	28"/30"	28"/30"	26"/28"	26"/28"	26"	24"	24"	24"	24"	24"
<b>2. Approx. Seat Length</b>				← 16" →				← 18" →			← 20" →	

Patient under 5'2" Hemi (Lowered) height, 5'4" with 2" cushion  
Patient in-between size, go down for men and up for women

**This chart is approximate, please measure for EXACT fit!**