

Trend QuickForm Hip

DCT-1686 The Trend QUICKFORM hip brace is designed to provide exceptional support and comfort for postoperative hip rehabilitation. It features a low-profile belt, optional hot/cold pack pocket, telescoping hip bar, adjustable range of motion dial, and an unloading hinge for precise recovery support. Developed by experts, it ensures advanced, patient-focused hip support.

Comfortable Design

The low-profile, moisture-wicking belt ensures a secure fit while keeping patients cool and comfortable during use. This design promotes long-term wear without compromising support.



Highlights

- **Therapeutic Integration:** Includes a pocket for optional hot/cold packs to aid in recovery.
- **Customizable Support:** Telescoping hip bar and range of motion dial allow for tailored adjustments.
- **Expert-Driven:** Designed by leading professionals for advanced patient compliance and care.

Indications

- Hip labral repair and resection
- Hip arthroscopy
- Hip Capsulotomy
- Femoroacetabular Impingement
- Hip support following surgery or injury

Recovery Optimization

The unloading hinge provides precise adjustments to offload pressure, helping to support the hip joint and promote faster recovery. It's ideal for meeting specific rehabilitation needs.

PDAC Approved: L1686

W trend-med.com P 800.441.2591 E info@trend-med.com

