



## VIVE OFFSET CANE MOB1002 HCPCS: E01003

The Offset Cane by Vive provides additional support and stability while walking by centering the user's weight over the cane. This allows for greater balance and support, while reducing hand and wrist fatigue. Fitted with a non-slip rubber cane tip, the offset cane can be used while navigating any surface. The cane height is also adjustable, for a comfortable and customized fit that can be used in either hand.

### WHAT IS INCLUDED

- 1x Offset cane

### MATERIALS:

- Anodized Aluminum
- Rubber Grip

### COLORS:



## DESIGNED FOR EASE OF USE AND PORTABILITY



# vive health

portal.vivehealth.com



**ADDS SUPPORT AND STABILITY**



**CUSTOMIZED SIZING**



**COMFORTABLE GRIP AND STRAP**

## **OFFSET HANDLE FOR EVEN WEIGHT DISTRIBUTION AND GREATER WRIST SUPPORT**

Ergonomically designed, the offset handle centers the individual's weight over the base of the cane for better balance and support. The offset hand grip also reduces hand fatigue for a comfortable experience.

## **ADDED STABILITY AND TRACTION ON ANY SURFACE**

Safely navigate over any surface with the slip-resistant rubber cane tip that is also non-marking. The offset cane is also compatible with most replacement cane tips.

## **COMFORTABLE ERGONOMIC GRIP**

Fitted with an ergonomically contoured hand grip, the cane handle creates a secure, yet comfortable grip perfect for daily use. The handle also includes a wrist strap for added convenience.

## **CORROSION-RESISTANT LIGHTWEIGHT FRAME**

Constructed with a durable anodized aluminum, the offset cane is exceptionally lightweight and corrosion resistant. Dual security features include a push button adjustment pin with a locking ring to ensure that the cane stays at the set height without slipping or rattling. The versatile cane is adjustable from 29" to a maximum of 38" and supports up to 250 pounds.

## **VIVE GUARANTEE**

Unconditional lifetime guarantee so you can purchase now with confidence.

**vive health**

8955 Fontana Del Sol Way  
Naples, FL 34109  
1-800-487-3808  
www.vivehealth.com

SR23- 1R1



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com