

vive
vivehealth.com

Distributed by

vive
health

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com

ANKLE BRACE

SUP1038



To see all the FAQ's,
visit vhealth.link/rx0

vive
vivehealth.com

QUICK START GUIDE

031023

HOW TO USE

Step 1:

Insert your foot into the brace so the Vive logo is over the front of your shin.



Step 2:

Fasten the top straps and adjust to your comfort level.



Step 3:

Attach provided ankle extender straps if needed. Extenders can be cut to desired length for a custom fit.



Step 4:

Wrap the bottom straps around your ankle and fasten them securely.

